



- non-invasive flow measurement of liquids
- for use on flexible tubing
- available for all common tube sizes
- hanging installation
- customer-specific calibration (seven calibration tables per sensor)

## Technical Specification

<b>Size (H x W x D) &amp; Weight (depending on ID)</b>	25 x 33 x 45 mm; 137 g 27 x 38 x 51 mm; 153 g
<b>Housing and Lid Material</b>	epoxy resin, aluminum, brass
<b>Cable Length</b>	2.9 m ±5 cm
<b>Safety class (IEC60601-1)</b>	BF (only in combination with the SonoTT™ Ultrasonic FlowComputer)
<b>IP-Code</b>	sensor head and cable: IP67 sensor connector: IP67 in mated condition
<b>Compatible Flow Meters</b>	SonoTT™ Ultrasonic FlowComputer
<b>Connector plug</b>	D-Sub
<b>Expected Product Life</b>	10 years

## Accuracy\* and Drift (in Combination with a SonoTT™ Ultrasonic FlowComputer)

Flow Range	Accuracy
0 < flow < Qmin	not defined
Qmin to 1000 ml/min	±100 ml/min
> 1000 ml/min	±7 % of the value ±30 ml/min

### \*Please note:

- The accuracy stated above can only be guaranteed if the parameters throughout the measurement are the same as the ones the sensor was adjusted and calibrated for (e.g. medium type, medium temperature, tube size, tube material).
- Qmin refers to the minimum flow value for which the accuracy is specified.
- Qmax refers to the maximum flow value for which the accuracy is specified.

## Adjustment and Calibration

<b>Recommended tube type</b>	flexible, non-reinforced tubing, e.g. silicone, PVC
<b>Medium type</b>	liquids such as water, purified water, salines solutions, blood, blood substitutes, and most dialysis, infusion, and irrigation solutions <b>Please note:</b> Due to safety and hygienic reasons, the adjustment is carried out using water.
<b>Medium Operating Temperature</b>	4 °C to 41 °C (40 °F to 105 °F)
<b>Calibration Tables</b>	up to seven calibration tables can be stored to each sensor plug

# SonoTT™ Clamp-On Transducer

## Range of SonoTT™ Clamp-On Transducers and Their Flow Measurement Range

Type	Qmin**	Qmax**	Tube Size (inner diameter (ID) x wall thickness (WT))
CT 6.8 mm	100 ml/min	±6 000 ml/min	6.8 mm outer diameter
CT 3/16" x 1/16"	100 ml/min	±6 000 ml/min	3/16" x 1/16"
CT 1/4" x 3/32"	120 ml/min	±8 000 ml/min	1/4" x 3/32"
CT 1/4" x 1/16"	120 ml/min	±8 000 ml/min	1/4" x 1/16"
CT 3/8" x 1/16"	150 ml/min	±10 000 ml/min	3/8" x 1/16"
CT 3/8" x 3/32"	150 ml/min	±10 000 ml/min	3/8" x 3/32"
CT 1/2" x 3/32"	300 ml/min	±20 000 ml/min	1/2" x 3/32"

### \*\*Please note:

While a measurement is generally possible for the flow range specified above, the accuracy is only defined for the flow values stated under "Accuracy in Combination with the SonoTT™ Ultrasonic FlowComputer".

## Ambient Conditions during Transport, Storage, and Operation

### Transport and Storage

Atmospheric Pressure	70 kPa to 106 kPa
Temperature Range	-20 °C to 55 °C (-4 °F to 131 °F)
Relative Humidity	10 % to 96 % (non-condensing)

### Operation

Atmospheric Pressure	70 kPa to 106 kPa
Operating Altitude	up to 3 000 m (9 842 feet)
Temperature Range	10 °C to 40 °C (50 °F to 104 °F)
Relative Humidity	10 % to 96 % (non-condensing)

**em-tec**

em-tec GmbH  
Lerchenberg 20  
86923 Finning, Germany  
P: +49 8806 9236 0  
F: +49 8806 9236 50  
E: em-tec-info@psgdover.com

em-tec.de

D 162-703 SonoTT™ Clamp-On Transducer- Technical Data Sheet - V5.0 | EMT-20002-T-01-A4  
© 2022 PSG, a Dover company

em-tec reserves the right to modify the information and illustrations contained in this document without prior notice.  
This is a non-contractual document.



Where Innovation Flows