
Capacitive Level Sensor (CLS) & Level Sensor Pad (LSP)

Copyright

All rights reserved, especially with respect to reproduction, distribution and translation. These Operating Instructions or any part thereof must not be reproduced, saved, processed, duplicated or distributed without the written consent of em-tec GmbH.

© Copyright em-tec GmbH Finning

Subject to technical changes

Owing to our policy of continuous product development, the illustrations and technical data contained in this document may differ slightly from the current version of the device.

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Table of Contents

1	General	23
1.1	Information about these operating instructions.....	23
1.2	Environmental protection.....	23
1.3	Notice concerning compulsory registration.....	23
2	Safety	24
2.1	Intendend use.....	24
2.2	Compatibility restrictions.....	25
2.3	General safety information for clinical use.....	25
2.4	Symbols.....	27
3	System description	28
4	Integrating into a base unit	29
4.1	PIN assignment.....	29
4.2	Signal output characteristics.....	30
4.3	Test procedure for checking functionality.....	30
4.4	Risk management for CLS/LSP device combination for a basic device (according to EN ISO 14971, IEC 60601-1, IEC 60601-1-2).....	31
4.5	Electrical and mechanical safety.....	32
4.6	Electromagnetic compatibility.....	32
5	Using CLS and LSP	33
5.1	Preparation.....	33
5.2	Setup.....	34
6	Cleaning and disinfection	36
7	Technical specifications of the CLS and LSP	37
7.1	Capacitive Level Sensor.....	37
7.2	Level Sensor Pad.....	37
8	Reorders	38
9	Support und Service	38

1 General

1.1 Information about these operating instructions

These operating instructions are for you to familiarise yourself with the features of the Capacitive Level Sensor (CLS) and the Level Sensor Pads (LSP) that go with it. Read these operating instructions carefully before use when integrating the CLS and LSP into your electrical device or medical product. The information for clinical users can serve as the basis for operating instructions on the clinical use of the CLS/LSP as an accessory.

Definitions

Indicates general imminent danger that may result in a person's death or severe injury.

Indicates potential danger that may result in a person's injury.

1.2 Environmental protection

Packaging

Packaging materials are made from environmentally friendly materials. The packaging materials will be disposed of by em-tec upon request.

Disposal

For the best utilisation of the raw materials, the project and its components and accessories should not be disposed of with household waste. All parts must be collected separately from household waste and disposed of in an environmentally responsible way in accordance with local regulations.

- Before disposal, decontaminate all parts according to the procedure applicable in the clinic.
- If you have questions about disposal, please contact em-tec GmbH's service department at: **service@em-tec.de**
- It may only be brought to the appropriate recycling facility if there is no risk of potential infection from electrical and electronic waste.

1.3 Notice concerning compulsory registration

Any major incidents connected to the product in any way must be immediately reported to the manufacturer and the appropriate authority of the member state where the user and/or the patient reside in.

2 Safety

2.1 Intended use

Proper intended use of the CLS/LSP

The CLS is intended as an electronic sensor component to integrate into electrical or electronic devices and medical devices, referred to hereafter as “base units”. The base unit must have a corresponding electronics system that supplies the sensor with the specified DC voltage through the connector and records both the current consumption and signal voltage at the interface and processes them accordingly. The reusable CLS is used in conjunction with the single-use LSP for monitoring fluid levels, such as blood or saline solutions, in a hardshell or bag reservoir. The CLS contains the sensor electronic system and emits changes in the electrical signal when a defined value does not reach the fill level. The monitored area is determined by the position of the LSP.

Intended users of the integrated CLS/LSP

For reasons of patient safety, the product may only be used by trained medical professionals, like a perfusionist or a trained nurse.

Intended environment of the integrated CLS/LSP

The product is intended for use in a clinical environment, away from strong magnetic and electromagnetic interfering fields (high-frequency surgery, magnetic resonance therapy, mobile phones).

Intended patients of the integrated CLS/LSP

The product can be used with all patients regardless of age, body weight and sex. The product is not for direct use on patients, and should not come into direct or indirect contact with the patient.

Intended medical use of the integrated CLS/LSP

The CLS was designed with the LSP as an independent accessory part for monitoring fluid levels (standard saline solutions or other aqueous solutions, and blood) in a reservoir as part of extracorporeal circulation. The safety and performance as an accessory were tested and proven with representative heart-lung machines and support systems. It can generally be used for solid polycarbonate, flexible polyvinyl chloride (PVC) or glass reservoirs. If the fill level does not reach a defined limit, a defined electrical signal change will be emitted.

Contraindications

The device must be used in the scope described in these operating instructions and for its intended purpose. Furthermore, the device may not be used for the following purposes:

The device is not intended for home care.

Additional contraindications and application restrictions can be found in the risk management for the full system, which should be integrated into the CLS/LSP.

2.2 Compatibility restrictions

Due to an increase in standardised requirements, the electronic system of the CLS has been updated to meet the current state of the art, so it is not intended to be integrated into base units that only match the signal characteristics and electrical properties of older versions of the CLS.

The company that distributes the base unit into which the CLS/LSP is integrated must refer the clinical user to the proper identification of the compatible and/or incompatible CLS variant(s) for the base unit, because due to the mechanical compatibility of the various connectors, it is possible that the device could be paired with the incorrect base unit.

The use of extension cords or adapters is only permitted if they have been tested and specified by em-tec and the manufacturer of the base unit for this purpose.

EN

2.3 General safety information for clinical use.

- The information for clinical use serves to minimise the risks associated with extracorporeal circulation and are based on *ISO/TS 23810* "Cardiovascular implants and artificial organs – Checklists for use of extracorporeal circulation equipment", risk management and market observation for the CLS/LSP, among other sources.
- Before using the system, read these operating instructions and the operating instructions for all approved single-use products and accessories that are used. In particular, take note of the operating instructions for the base unit and the container to be monitored regarding positioning and use of level monitoring.
- Use the system and its accessories in accordance with the operating instructions and proven medical practice. Also take note of the relevant current guidelines and handling recommendations of medical expert associations and organisations.
- Install the level monitoring system according to the instructions for the base unit.
- Only start the application with fully functional level monitoring. Before use, perform a function test according to the instructions for the base unit.
- The adhesive strip on the LSP adheres very tightly to clean surfaces; residue will remain when it is removed. The LSP cannot be reused, as its adhesive strength will be lost. A replacement LSP should not be stuck onto the remnants of an adhesive strip, as this can cause malfunctions.
- Repeat the function test if the LSP is replaced or its position is changed. If the function is defective, use a new LSP or replace the CLS.
- Keep a tested and functional CLS and LSP on hand so that in case of a breakdown, damage to the level monitoring system or a change in the system, the level monitoring system can be recovered.
- When handling blood and other potentially infectious fluids, observe the safety regulations and hygiene protocol in place in the clinic.
- Observe the requirements for position of use, operating position and positioning of the heart-lung machine or its support system.

- Plugging in and unplugging wet sensor connectors can cause measured values to become corrupted. In some circumstances, the removal of a sensor connector from the heart-lung machine or its support system will not be recognised and no alarm will be triggered. Wet sensor connectors and safety caps must be dry before plugging in or unplugging.
- A short circuit in the CLS cable can cause malfunctions. There is the risk that a malfunction in the level monitoring system and failure to meet the level limit will not be detected.
- The level monitoring system with connected CLS may only be operated with an LSP when the LSP is fully adhered to the container to be monitored. A CLS with connected LSP that is not affixed to the container to be monitored may detect falsely metallic objects or surfaces (e.g. RFID labels, CLS cables or a second CLS or LSP) or human bodies within range of the LSP as fluids.
- Keep components safe from moisture.
- Protect the CLS from mechanical stress and damage. Do not hit it.
- Do not let the CLS hang loose. If you do not need the CLS, mount it securely. One way to do this, for example, is to stick an LSP in a secure spot where the CLS can be mounted.
- Do not let any hoses or cables hang loose.
- Do not let components fall. Components that have been exposed to mechanical shocks (such as from falling) must be taken out of service and inspected by the authorised technical service.
- The level sensor must be cleaned and disinfected after each use and visually inspected for visible damage.
- Accuracy may vary if CLS/LSP is attached to a moving surface.
- CLS/LSP must not be used for controlling, monitoring or directly influencing the performance of active implantable devices.
- Absence or degradation of level sensing and its components CLS and LSP will NOT result in an unacceptable risk for the patient, user or environment; (SELV, uncritical material, class B or better).

- Make sure that the cables for the CLS are not excessively bent, as this can damage the ferrite core.
- The device may not be used with other fluids that are not water-based, such as hydrocarbons, oil etc.
- CLS-LSP may not be used on a reservoir that is smaller than the LSP.
- The CLS cannot be combined with level pads other than those approved by em-tec.
- Hardshell tanks or containers whose geometry does not permit the LSP to be fully horizontally adhered or on which the adhered LSP does not permit fill level to be monitored visually may not be used.

2.4 Symbols

Symbol	Meaning
	Warning
	Attention
	Observe operating instructions
IPX 1	Protection Class under <i>IEC 60529</i> : protection against drops of water falling vertically.
	Separate collection of electrical and electronic devices in accordance with Directive <i>2012/19/EU</i> : The product may not be disposed of with household waste. It must be collected separately from household waste and disposed of in an environmentally responsible way in accordance with local regulations.
	Manufacturer according to Medical Device Regulation (<i>EU</i>) <i>2017/745</i> .
	Date of manufacture: Date on which the product was manufactured (year).
	The product fulfills the general requirements of Medical Device Regulation (<i>EU</i>) <i>2017/745</i> . The digit 0123 is the code number of the responsible notified body TÜV Süd Product Service GmbH, Munich
	Serial number
	Order number
	Temperature limit
	Humidity, limit
	Air pressure, limit
	Fragile!
	Keep dry

Symbol	Meaning
	Medical Device
	MR Unsafe! Item should not enter the MRI scanner room. Patients with MR Unsafe devices should not be scanned.
	Active mobile phones prohibited

3 System description


Fig 1: CLS and LSP – display of CLS clamping direction.

The level monitoring system consists of the ① “Capacitive Level Sensor” (CLS) and the ② single-use product “Level Sensor Pad” (LSP).

The CLS, combined with the LSP, is an independent accessory part for monitoring fluid levels (standard saline solutions or other aqueous solutions, and blood) in a reservoir. It can generally be used for solid polycarbonate, flexible polyvinyl chloride (PVC) or glass reservoirs.

The CLS is clamped onto the end of the LSP and fixed in place by the force of a built-in clip. The contact pins are electrically connected to the capacitive aerial on the LSP.


Fig. 2: Example of how to position the LSP on a reservoir.

If the horizontal level of fluid within the area ③ determined by the LSP’s position changes, the signal level status at the output changes, which is transmitted to a base unit.

The sensors have not been disinfected and are not sterile. The LSPs are only delivered in nonsterile condition and should be kept in the shipping bag until they are used, as this is where the information on the batch number, expiry date and, if applicable, UDI are found and the LSP is protected from environmental factors.

Classification under MDR (EU Directive 2017/745)

The CLS medical device with its LSP accessory is assigned to MDR Class IIb in accordance with Annex VIII, using Rule 9 and Rule 1. The CLS is an active device used to monitor the performance of active therapeutic devices (eg, heart-lung-machines). The CE certification from TÜV SÜD (MDR) for em-tec GmbH is based on this classification.

A copy of the current version of the CE declaration of conformity, including RoHS, and a copy of the CE certificate are not included in this operating manual, but are available from em-tec GmbH by request.

EN

4 Integrating into a base unit

4.1 PIN assignment

As explained in section 3 System description, the system consists of one CLS, which contains the electronic system, and the LSP. The CLS is clamped onto the end of the LSP. The contact pins are connected to the capacitive aerial on the LSP (fig. 2, section 3)


Fig. 3: PIN assignment for the CLS.

The cable is shielded. The shield is electrically connected to the connector housing and the CLS housing so that the potential on the connector is also in contact with the anodised CLS housing. Electrical measurements, such as contact voltage, should be measured on the cable inlet, as the anodised coating forms an insulating barrier.

- There are other versions of the CLS on the market that differ significantly in their electrical layout, electrical characteristics, signal behaviour and/or electromagnetic compatibility. However, these variants use the same type of connector with the same PIN assignment, amount of contacts and the same size. Make absolutely sure that you are using the right CLS variant for your application or to integrate into your base unit. The different versions can be distinguished from one another by means of their ID (item number) or the operating instructions provided with each device. Mixing up versions of the CLS will not cause any danger to patients or users.
- Some combinations of connector and socket from different manufacturers (Fischer and Lemo) are mechanically compatible. The CLS will not work if the combination is incompatible.
- The threshold voltage depends on the supply voltage.

4.2 Signal output characteristics

The signal is emitted as an analog voltage signal. When the fill level is above the limit set, a signal of 0 volts will be emitted. If the fill level is within the limit range ③ or underneath, a voltage of 3.5 volts \pm 5% will be emitted (fig. 4).


Fig. 4: Example of how to position the LSP on a reservoir.

4.3 Test procedure for checking functionality

Carefully test the system after it has been integrated into the base unit. During starting and normal operation, the function and the connection of the CLS and LSP should be monitored by the base unit so that malfunctions, electrical short circuits and detachment can be detected, e.g. by measuring the current consumption or analog output.

em-tec recommends the following tests to verify the device's integration into the base unit. It is required that these tests be performed before completing integration and releasing the CLS/LSP for the base unit:

CLS and LSP fluid level detection test	
1	Use an empty reservoir or a reservoir with a fluid level below the LSP.
2	When the CLS is connected to the LSP, the base unit must respond and display the correct sensor status as intended.
3	Increase the level of fluid in the container so that the fill level is between the lower and upper limits of the alarm activation range.
4	The base unit must respond and display the correct sensor and fluid status as intended.
5	Increase the level of fluid in the container so that it is above the LSP.
6	The base unit must respond and display the correct sensor status with full fill level as intended.
7	Reduce the level of fluid in the container again so that it is below the LSP.
8	The base unit must respond and display the correct sensor status for the low level of fluid as intended.

CLS and LSP connection test	
1	Fill the container until the level of fluid is above the LSP.
2	When the CLS is connected to the LSP, the base unit must respond as designed and display the correct sensor status and fluid level as intended.
3	When the CLS is not connected to the LSP or the electrical connection to the LSP is faulty, the base unit must respond and show that fill level monitoring is not active, as intended.

4.4 Risk management for CLS/LSP device combination for a basic device (according to EN ISO 14971, IEC 60601-1, IEC 60601-1-2)

EN

The electronic system for the CLS is based purely on electronic components, without the use of software. Some error statuses are, as described, signalled by changing the current consumption or the voltage of the signal output.

The company that distributes the combination of CLS/LSP and base unit is responsible for risk management under the environmental and operation conditions provided to them, as these cannot be accounted for by em-tec GmbH when designing, developing and producing the CLS/LSP. Monitoring and adequate response to the specified electrical characteristics of the CLS/LSP (change of signal voltage, power consumption) must be monitored by the base unit. Known risks and risk reduction measures resulting from the combination and interaction of CLS/LSP that need to be assessed can include the following:

- Breakdown or change in sensitivity of level detection and its impact on safety and/or essential features of the base unit.
- Function test to detect damaged or soiled products or breakdowns in functionality of the level detection system, incorrect positioning and faulty behaviour in the base unit before a patient is hooked up to extracorporeal circulation.
- Provision of replacement parts for the CLS and LSP.
- Usability according to *IEC 60601-6* or *IEC 62366-1*.
- If two bags are attached one after the other, make sure that the CLS/LSP only detects the limit range for the bag in front for which the fill level should be monitored. In this instance, the integrator and, if necessary, the user will need to test this situation.
- User information, product labelling and correct product identification/selection.
- The manufacturing process for the CLS/LSP combination with the base unit must ensure that the specified electrical properties are guaranteed at the interfaces and that the system responds as designed.
- The minimum/maximum permissible limits must be set for reservoirs in order to prevent the risk of the reservoir running over or running empty as is not permitted, while taking into account the system characteristics and the volumes of fluid that are supplied or discharged. Varying positions or fluid surfaces that fluctuate or aren't horizontal must be taken into account for this.
- Interference due to the use of two or more CLS units in close proximity to each other or other products.
- The manufacturing process for the base unit ensures that the required electrical power is guaranteed. Quality management *ISO 13485* is a prerequisite.

4.5 Electrical and mechanical safety

The intended transport, storage and application conditions in addition to the interface requirements for integration into a base unit are specified in this document.

Modifications and repairs are not permitted. The CLS has no components that must be maintained or replaced. The sensor is operated with safety extra-low voltage (9V DC). The company that distributes the CLS/LSP combination with the base unit must ensure that the standardised requirements in the *IEC 60601* and *ISO 80801* standards concerning electrical shock (MOP) for users (MOOP) and patients (MOPP) and touch and leakage currents (at least Class B) are complied with. The unit is not designed to come into contact with the patient and this should be avoided. Please note that the black anodised coating has additional insulation, while the potential on the cable inlet is the same as that on the connector due to the shield.

The sensor has a single channel and provides the electrical characteristics described, so it can detect potential error statuses or conditions in which level monitoring is not possible.

The voltage supply and signal evaluation of the base unit must be safeguarded so that the first instance of malfunction, such as a short circuit or broken line, does not cause unacceptable risk.

The electrical safety and performance was metrologically measured and confirmed for the CLS and LSP as accessories in combination with representative heart-lung machines and support systems in an accredited testing lab in accordance with *IEC 60601-1*.

The maintenance and necessary safety and function tests must be specified by the manufacturer of the base unit in accordance with *IEC 60601-1* and *IEC 62353* in order to guarantee the safety of the CLS during its service life. Tests may include visual inspection, the measurement of leakage currents and the consistency of the shield on the sensor housing. If tests such as insulations measurements are done incorrectly, however, this may result in damage to the sensor electronic system.

The tests are to be coordinated with em-tec.

Testing for mechanical stress is based on the base unit's classification. The system integrator must ensure that the CLS is not damaged when used as intended as per the base unit specifications.

The minimal distance between capacitive level sensors (when applicable) or to components of the base device have to be specified and tested for the device (e.g. heat exchanger).

4.6 Electromagnetic compatibility

The test from *IEC 60601-1-2*, 4th edition, can only be performed in conjunction with the base unit while level monitoring is active for it.

In order to meet the standardised requirements for electromagnetic compatibility with a base unit, the CLS cable can be equipped with optional qualified ferrite cores to act as measure against interference.

A function test is required before using the device on patients in order to detect damage on the electrical assembly in the CLS or the level monitoring electronics system in the base unit that may result from electrostatic discharge, or to identify environmental conditions that impact level monitoring, such as electromagnetic interference fields or conductive objects, especially in the area around the LSP's integrated aerial.

The immunity of the level detection to foreseeable electromagnetic interference fields such as those described in *IEC 60601-2-1*, for example, can only be assessed using the CLS in combination with the base unit.


In order to guarantee proper level detection, there must be sufficient distance from electromagnetic interference fields (such as from mobile phones) in the environment in which the device is used. When using medical applications with interference fields such as high-frequency surgery, it may be necessary to disable level monitoring and use other methods, such as visual monitoring.

The electromagnetic compatibility was metrologically measured and confirmed for the CLS and LSP as accessories in combination with representative heart-lung machines and support systems in an accredited testing lab in accordance with *IEC 60601-1-2*.

The implementation of public mobile radio standard 5G or the use of other wireless technologies, such as RFID, will entail new frequency ranges and field strengths that result in requirements for immunity that are not specified by the 4th edition of *IEC 60601-1-2*, although they must also be factored into the risk management conducted by the company that distributes the base unit.

5 Using CLS and LSP

5.1 Preparation

Lay the LSP on a flat surface and peel the protective film away from the adhesive side of the sensor pad. Only use LSPs that are not expired or damaged. Otherwise, the sticker on the LSP may not be sticky enough, which can result in malfunctions or false alarms. Avoid contact with the bond area of the Level Sensor Pad. The adhesive tape will be damaged if the LSP is removed from the reservoir.

Fig. 5: Remove the protective film on the LSP

- Avoid flexing or bending the LSP too much. This can cause damage to the aerial and the alarm function may not work.
- The LSP's adhesive surface area measures 72 x 26 mm \pm 5%, and must be fully horizontally applied to the surface underneath it in order to guarantee the sensor characteristics as intended.
- In rare cases, the sensor may signal an empty bag when the fill level is above the LSP.
- Do not let cables, liquid lines (hoses) or liquids enter the area between the LSP nipples.
- The diameter of the reservoir must not be less than 80 mm.
- The wall thickness must not be more than 5 mm.


Fig. 6: LSP trigger range

The LSP's position defines the monitored fill level area. The fill level sensor generates a change in the level of the electrical output signal when the level of fluid moves within the triggering range ① of the LSP area. The medical base unit can then generate an alarm and/or trigger an action. It must be positioned horizontally with the tip pointing to the right so that the CLS cable hangs down vertically.

5.2 Setup


Fig. 7: Sample positioning for LSP trigger range

Fasten the LSP to the fluid container to be monitored that corresponds to the fill level at which the system should respond when the fluid level is not high enough (limit range ②, fig. 7). Also follow the manufacturer specifications and/or labelling on the reservoir. The LSP may only be attached and used in the filling area specified and permitted for the respective container.

The LSP should be attached in such a way that the CLS cable is not bent. The lower threshold of the alarm activation range must be above the critical fluid level.

CLS and LSP should be inspected for residual dirt, fluids or other contaminants, especially in the area of contact.

Test before use

Before using the capacitive level sensor with the CLS/LSP, an appropriate function test must be conducted by the user in conjunction with the base unit. Possible test steps for checking the sensor's functionality can be found in section 4.3, "Test procedure for checking functionality".

Removing/shutting off

Remove the CLS from the LSP after use. Dispose of the LSP with the components to which it is stuck. The CLS should be cleaned and disinfected after use. Do not store the device until it is completely dry. Avoid bundling or bending the cable too much when winding it up, especially around the tension relief device or the ferrite core.

Using multiple CLS/LSPs simultaneously

The user must ensure that the correct CLS/LSP is matched with the proper connection to the heart-lung machines or support system.

Mishandling/electrical behaviour

Problem	Meaning / solution
The current consumption exceeds the specified range	Contacts are shortened/bent and/or cable is faulty ⇒ Replace the CLS
No current consumption/ no analog signal	CLS is not connected ⇒ Reconnect the CLS or replace it Damaged and/or bent contacts ⇒ Replace the CLS
Analog output > 6.0 V	Short circuit between supply line and signal line ⇒ Replace the CLS
The analog output does not respond or does not respond adequately when the fill level is low	Inspect the CLS and LSP for wetness from fluids, moisture, overflow and dirt, especially around the contacts. ⇒ Clean the CLS and let it dry Liquid film, foam or residue inside the reservoir near the fill level sensor is interpreted incorrectly. ⇒ Place the new LSP in a better location
When fill level is high, the analog output does not respond or does not respond adequately	The reservoir used is incompatible with CLS and/or LSP. ⇒ Replace the container with one of a compatible type
The analog output does not respond to changes in fill level	Check whether the aerial contacts on the LSP are loose ⇒ Replace the LSP
The analog output level changes accidentally	Check whether there are cables or hoses from infusions - or extracorporeal hose systems touching or near the back of the LSP. ⇒ Eliminate source of distortion, ensure adequate distance Check whether residual fluids or residues are on the back of the LSP. ⇒ Eliminate cause of signal distortion ⇒ Replace the LSP with a new one if necessary
LSP is not recognised by the CLS	Check whether the aerial contacts on the LSP are loose ⇒ Replace the LSP

6 Cleaning and disinfection

- Make sure that cleaning agents and disinfectants do not enter the device or its accessories. Do not immerse the CLS in or spray it with liquids. Only use a damp, but not wet, cloth to clean and wipe off the CLS. Therefore, use the CLS/LSP system only outside the sterile area.
- It is not meant to be sterilised.
- Only use the alcoholic disinfectants specified and follow the manufacturer's instructions regarding application and dosage.

For cleaning the CLS, we recommend the use of common mild cleaning agents that are suitable for sensitive materials. To prevent layers of dirt from forming, residual dirt should be removed after each use. Use a lint-free cloth to clean the surface of the CLS.

Do not adhere the LSPs to freshly disinfected or cleaned surfaces; wait until the surface is dry.

Disinfection (or sterilisation) methods other than those specified must be specified and validated by the responsible manufacturer of the base unit, especially those concerning safety, functional and material compatibility or effectiveness.

The following alcohol-based disinfectant is recommended for the CLS:

Name	Manufacturer	Contact
Bacillol	Bode Chemie	Scientific Affairs E-Mail: KundenService-SiDa@bode-chemie.de

7 Technical specifications of the CLS and LSP

7.1 Capacitive Level Sensor

Housing dimensions	42 x 30.6 x 23 mm ± 5%
Housing material	Aluminium alloy, black anodised
Weight	28g ± 5g
Cable	Insulated PVC, 3 wires, shielded Length: 1970mm ± 30mm
Protective class	IPX1 (dripping water): according to EN60529, the electronic equipment of the CLS is protected from water penetration.
Connector type on cable	Lemo FFG.0B-303
Compatible socket type	Lemo 0B303
Voltage supply	DC 9V ±5%
Max. energy consumption	135 mW
Current consumption connected to LSP	10-13 mA
Current consumption not connected to LSP	<1.5 mA
Threshold voltage when fluid level is low	+3.5V ±5% DC
Wall thickness of the reservoir	max. 5 mm
min. diameter of the reservoir	min. 80 mm

Environmental conditions

Operating temperature	10°C to 40°C (50°F to 104°F)
Transport temperature	-20°C to 55°C (- 4°F to 131°F)
CLS storage temperature	-20°C to 45°C (- 4°F to 113°F)
Relative operating humidity	30-75%, non-condensing
Relative storage and transport humidity	10-96%, non-condensing
Operating altitude	Up to 3000m (9842 feet)
Atmospheric pressure during transport, storage and operating	70 to 106kPa

7.2 Level Sensor Pad

Dimensions	72 x 26 x 13 mm ± 5 %
Weight	3.5 ± 0.5 g
Plastic clip	Polyethylene with low density
Aerial	Aluminium foil
Adhesive strip on sensor strip	0.8mm large white, interlinked, closed-cell polyethylene foam, coated with rubber-based adhesive (latex-free)

Environmental conditions (difference from CLS)

LSP storage temperature	15°C to 25°C (59°F to 77°F)
LSP storage	in bag, do not remove until use, see LSP label for storage conditions
Other	All other specifications same as CLS

8 Reorders

Capacitive Level Sensor

- 1 CLS (9 Volt) Lemo Connector, includes operating instructions

Order number: M157-002

Level Sensor Pad

- Package with 100 LSPs and short guide (use-by date 2 years from production)

Order number: 12747

Using the contact information below, you can also learn about our ordering and shipping methods for yourself. Our product experts are happy to advise you on the possible uses and performance data of our products:

em-tec GmbH
Lerchenberg 20
86923 Finning
GERMANY

Phone: +49 8806 9236 0
Fax: +49 8806 9236 50
E-Mail: info@em-tec.de
Internet: www.em-tec.de

9 Support und Service

If the recommended measures listed in the “Troubleshooting” section are not successful, please contact our customer service. Take down the serial number of the CLS and/or the batch number of the LSP or, if available, the UDI, before contacting the department. Please provide as much information about the combination of devices (e.g. base unit, container type, fluid) and the conditions under which they are used as possible.

Please call or e-mail to reach our service department:

Phone: +49 8806 9236 0

E-Mail: service@em-tec.de

If a device needs to be sent in for repair, please request an RMA (Return Material Authorisation) number from the service department and enclose a detailed description of the error and/or problem as well as name and number/e-mail address of contact in the shipment. You can find the form at the end of this manual.

Please use the original packaging for the shipment, as it is specially designed for transporting this device.

Shipping address:

em-tec GmbH
Service Department
Lerchenberg 20
86923 Finning
GERMANY

RETURN MATERIAL AUTHORISATION FORM

		Serial No.:	
--	--	--------------------	--

Customer:			
Company:			Contact:
Address:			Phone:
City:			Fax:
State:	Zip:		E-mail:

Delivery address for return shipment to customer: same as above: • if not specify below:

Company:			Contact:
Address:			Phone:
City:			Fax:
State:	Zip:		E-mail:

Reason for the RMA request (description of the malfunction in detail):

Was there an injury or death? Yes No Unknown

Could this defect lead to an injury or death? Yes No Unknown

Name: _____ **Date:** _____ **Signature:** _____

Instructions:

- Please complete this form, print and sign it and attach it to the shipment.
- Please send the shipment to the em-tec Service department:
em-tec GmbH
 Service department
 86923 Lerchenberg
 20 Germany

Contact info: phone +49 8806 9236 0 / Fax +49 8806 9236 50 / e-mail: service@em-tec.de

- Declare on the delivery note that the unit will only be imported temporarily for repair (important for customs clearance of shipments arriving from outside the EU).
- Please also note the serial number and the RMA number on the delivery note.
- Before packaging, please make sure that the goods are cleaned and disinfected.
- Make also sure that the goods are carefully packed, if possible in the original package.
- Transport insurances with sufficient coverage can be helpful to cover transit damages.

Special Note: If a price quote (±20%) shall be provided before repair, please check the box. Please note that a handling fee of 100.00 EUR needs to be charged for the quote in addition to costs for freight and insurance (terms of delivery: EXW).	<input type="checkbox"/> Provide quote before repair
--	---

EN