
# SonoTT™ Clamp-On Transducer

## Operating Instructions

## **Copyright**

All rights reserved, especially with respect to reproduction, distribution and translation.

These operating instructions or any part thereof must not be reproduced, saved, processed, duplicated or distributed without the written consent of em-tec GmbH.

© Copyright 2022 em-tec GmbH Finning. -All rights reserved.

## **Subject to Technical Changes**

Owing to our policy of continuous product development, the illustrations and technical data contained in this document may differ slightly from the current version of the device.

## **Legal Manufacturer**

em-tec GmbH  
Lerchenberg 20  
86923 Finning  
GERMANY

Phone: +49 8806 9236-0  
Fax: +49 8806 9236-50  
E-mail: [em-tec-info@psgdover.com](mailto:em-tec-info@psgdover.com)  
Internet: [www.em-tec.de](http://www.em-tec.de)

## Table of Contents

English.....	5
Deutsch.....	26

EN

DE

## Table of Contents

<b>1</b>	<b>System Description.....</b>	<b>5</b>
<b>2</b>	<b>Description of the Measurement Principle.....</b>	<b>6</b>
<b>3</b>	<b>General Safety Information, Symbols and Units .....</b>	<b>7</b>
3.1	Symbols Used in these Operating Instructions .....	7
3.2	Symbols on the SonoTT™ Clamp-On Transducer and on Packaging.....	7
<b>4</b>	<b>General Information Concerning these Operating Instructions.....</b>	<b>8</b>
<b>5</b>	<b>Intended Purpose, Restrictions and Limitations .....</b>	<b>8</b>
5.1	Intended Purpose.....	8
5.2	Contraindications.....	9
5.3	Usage Restrictions and Limitations.....	9
5.4	Liabilities and Responsibilities.....	10
5.5	Notice Concerning Compulsory Registration.....	10
5.6	Safety Instructions.....	10
5.7	Electromagnetic Compatibility .....	10
5.8	Ultrasonic Safety Regarding Acoustic Output.....	11
5.9	Biocompatibility .....	11
5.10	Maintenance and Service.....	12
<b>6</b>	<b>Installation, Use, and Disassembly of the SonoTT™ Clamp-On Transducer.....</b>	<b>13</b>
6.1	General Safety Information .....	13
6.2	Connecting the SonoTT™ Clamp-On Transducer to the SonoTT™ Ultrasonic FlowComputer.....	13
6.3	Attaching the SonoTT™ Clamp-On Transducer to the Tube.....	14
6.4	Use of the SonoTT™ Clamp-On Transducer.....	15
6.5	Disassembly and Storage of the SonoTT™ Clamp-On Transducer .....	16
6.6	Cleaning and Disinfection of the SonoTT™ Clamp-On Transducer .....	17
6.7	Changing the SonoTT™ Clamp-On Transducer .....	17
<b>7</b>	<b>Adjustment and Calibration of the SonoTT™ Clamp-On Transducer.....</b>	<b>18</b>
7.1	Adjustment.....	18
7.2	Calibration .....	18
<b>8</b>	<b>Troubleshooting .....</b>	<b>19</b>
<b>9</b>	<b>Environmental Protection and Disposal.....</b>	<b>22</b>
<b>10</b>	<b>Technical Specifications of the SonoTT™ Clamp-On Transducer.....</b>	<b>23</b>

## 1 System Description

The SonoTT™ Clamp-On Transducer, which is used in combination with, and powered by, the evaluation device SonoTT™ Ultrasonic FlowComputer, measures the volumetric flow rate of fluids based on the ultrasonic transit time method. The SonoTT™ Clamp-On Transducer encloses flexible tubes via a simple click-fastening and operates non-invasively without physical contact to the medium, allowing it to meet even strict hygienic requirements. In combination, the SonoTT™ Clamp-On Transducer and the SonoTT™ Ultrasonic FlowComputer represent an independent system.

Established fields of application on extracorporeal tubing circuits for the SonoTT™ Clamp-On Transducer are, for example:

- heart-lung machines (HLM)
- extracorporeal membrane oxygenation (ECMO)
- organ perfusion
- left and right ventricular bypasses
- cardiopulmonary support/bypasses
- venous extracorporeal support (e.g. for liver transplantation)
- shunts

Figure 1: SonoTT™ FlowMeasurement System

### Classification under MDR (EU Directive 2017/745)

The SonoTT™ Clamp-On Transducer is assigned to MDR Class I in accordance with Annex VIII, using Rule 13. The SonoTT™ Clamp-On Transducer is an active device used to monitor the flow in extracorporeal tubing circuits together with the SonoTT™ Ultrasonic FlowComputer. The Declaration of Conformity for the MDR is based on this classification.

While a copy of the current version of the CE Declaration of Conformity, including RoHS and REACH, and a copy of the CE certificate itself are not included in these operating instructions, they are available from em-tec GmbH upon request.

## 2 Description of the Measurement Principle

The SonoTT™ Clamp-On Transducer is based on an acoustic measurement principle and utilizes the transit time method to determine the flow. A typical transit time flow measurement system incorporates two piezo ceramics that act as both ultrasonic transmitter and receiver. The measurement system, consisting of the SonoTT™ Ultrasonic FlowComputer and the SonoTT™ Clamp-On Transducer, operates by alternately transmitting and receiving an ultrasonic pulse between the piezo ceramics and measuring the transit time difference that it takes for the pulse to travel between them.

This method determines the transit time difference of ultrasonic signals through a measuring section and is used to estimate volumetric flow rates. There are at least two piezo ceramics required in order to send and receive ultrasonic pulses both with and against the flow direction. For a higher accuracy, the SonoTT™ Clamp-On Transducers incorporate two pairs of piezo ceramics, i.e. four piezo ceramics in total, which are arranged in an X-configuration around the tube containing the measured flow.


Figure 2: Sensor Structure

When sending ultrasonic signals through the measuring section, the transit time difference depends on the flow direction of the medium:

- The ultrasonic sound signals that are sent along the flow direction and volume flow of the medium, i.e. downstream, need less time to travel the measuring section than
- the ultrasonic sound signals that are sent against the flow direction, i.e. upstream.

Consequently, the transit time difference is measured for every pulse. The difference between upstream and downstream measurements is proportional to the volumetric flow rate of the liquid in the tube. Therefore, by accurately measuring the difference between upstream and downstream transit time  $T_{up}$  and  $T_{down}$ , it is possible to estimate the flow velocity. Subsequently, the volumetric flow rate is computed from the inner cross section area of the tube and the instrument coefficient.

### 3 General Safety Information, Symbols and Units

#### 3.1 Symbols Used in these Operating Instructions

Symbol	Meaning
	Warning! This safety symbol precedes critical information that must be strictly observed in order to prevent injuries and fatal hazards. This warning symbol is the most important safety symbol.
	Caution! Important information regarding correct handling. Must be observed and strictly adhered to. If this information is not observed, malfunction or damage of the product and its surroundings may occur.

EN

#### 3.2 Symbols on the SonoTT™ Clamp-On Transducer and on Packaging

Symbol	Meaning
	Read these operating instructions carefully before use. Keep them in an easily accessible location for future reference. Replacement copies of these operating instructions are available from the manufacturer.
	Do not dispose of this device together with domestic waste! The device as a whole as well as any parts must be disposed of in accordance with WEEE Directive and national legislation.
	Manufacturer: em-tec GmbH · Lerchenberg 20 · 86923 Finning · Germany
	Date of manufacture
	Serial number
	Order number
	Temperature limit during storage and transport.
	Moisture limit during storage and transport.
	Air pressure limit
	Classification Applied Part: BF according to EN 60601-1
<b>IP67</b>	IP (=Ingress Protection) Classification

## 4 General Information Concerning these Operating Instructions

- Read these operating instructions carefully before starting up the system!  
→ This document is intended for the SonoTT™ Clamp-On Transducer and describes the use of the SonoTT™ Clamp-On Transducer in combination with the SonoTT™ Ultrasonic FlowComputer.
- The user (= the person who connects the SonoTT™ Clamp-On Transducer to the tubing circuit and measurement system) is responsible for any risk if the SonoTT™ Clamp-On Transducer is not connected correctly.
- Inspect your SonoTT™ Clamp-On Transducer for completeness and damage when unpacking it for the first time.
- The customer must ensure that the persons involved in the integration and use of the SonoTT™ Clamp-On Transducer and the system it is part of are adequately qualified. In addition, the information in these operating instructions must be followed.
- Do not drop the SonoTT™ Clamp-On Transducer. It might break or become damaged, which could impair the measurement.

These instructions contain important information concerning the safe and correct handling of the SonoTT™ Clamp-On Transducer. Read these operating instructions carefully before using the product and keep them in an easily accessible location. Familiarize yourself with all warning and safety information and observe them. It is the responsibility of the operator to ensure that the SonoTT™ Clamp-On Transducer is used, inspected and maintained in accordance with the operating instructions. This also applies for all subsequent revisions or instructions from the manufacturer. The manufacturer reserves the right to modify technical data without prior notice.

## 5 Intended Purpose, Restrictions and Limitations

### 5.1 Intended Purpose

The SonoTT™ Clamp-On Transducer (in combination with a SonoTT™ Ultrasonic FlowComputer) is indicated for the bi-directional volumetric measurement of the liquid flowing through tubing systems. The measurement principle is the ultrasound transit-time method.

The medical use of the device is appropriate for extracorporeal applications together with the SonoTT™ Ultrasonic FlowComputer such as:

- cardio-pulmonary bypass, hemodialysis, hemofiltration, plasmapheresis
- ECMO (ExtraCorporeal Membrane Oxygenation) applications
- perfusion, infusion, transfusion
- several shunt-applications

#### Intended Users of the SonoTT™ Clamp-On Transducer

For reasons of patient safety, the device may only be used by qualified and trained medical professionals.

### Intended Patients of the SonoTT™ Clamp-On Transducer

The product is intended for the following patient population:

- patients who are undergoing a surgery, supervised by medically trained personnel.
- patients who are connected to an extracorporeal circulation system such as an ECMO device.
- patients who are treated within an hospital environment, e.g. an intensive care unit or operating room.

The product is not intended for the direct use on patients and should not come into direct contact with them. In situations where skin contact cannot be ruled out or avoided, a pad, e.g. a gauze, should be placed between the skin and the sensor.

A malfunction of the device is causing no additional health risk to the patient.

### Intended Environment of the SonoTT™ Clamp-On Transducer

The SonoTT™ Clamp-On Transducer is designed for the continuous operation in intensive care units, operating rooms and similar medical environments.

### Clinical Benefits

Additional information for the attending physician on the flow behavior of liquids in a patient-bound tubing system.

The expected lifetime of the device is ten years.

## 5.2 Contraindications

This device and the accessories have been built and sold exclusively for the intended purpose mentioned above.

The SonoTT™ Clamp-On Transducer and the SonoTT™ Ultrasonic FlowComputer are explicitly **not intended**

- to allow direct diagnosis or monitoring of vital physiological processes.
- for the monitoring of vital physiological parameters, where the variation could result in immediate danger to the patient, for instance variations in cardiac performance, respiration, activity of the central nervous system.
- as sole base for decisions regarding therapeutic methods. The values issued by this device have to always be verified with other methods before any decisions regarding treatment are made.

It is the operator's responsibility to use the device as described in the operating instructions and in any of its following revisions. The use of clinical procedures and techniques lies in the physician's area of responsibility.

## 5.3 Usage Restrictions and Limitations

The SonoTT™ Clamp-On Transducer is explicitly **not intended**

- for the use outside of medical environments.
- for the use inside means of medical transport (e.g. ambulance car or helicopter).

## 5.4 Liabilities and Responsibilities

The user is responsible to use, check and maintain the SonoTT™ Clamp-On Transducer according to these operating instructions.

em-tec GmbH is neither liable nor responsible for any consequences arising from the use of the SonoTT™ Clamp-On Transducer that does not comply with the operating instructions or specifications in this document.

## 5.5 Notice Concerning Compulsory Registration

Any major incidents in any ways connected to the product must be immediately reported to the manufacturer and the appropriate authority of the member state in which the user and/or patient reside in.

## 5.6 Safety Instructions

If the SonoTT™ Clamp-On Transducer is part of a medical application, the user has to observe the following instructions:

- The measured values of the SonoTT™ Clamp-On Transducer, supplied through the corresponding evaluation device (i.e. the SonoTT™ Ultrasonic FlowComputer) only display the current flow rate (e.g. of an oxygenator or blood pump).
- These values are neither to be used for the evaluation of a patient's status nor for the basis of medical actions without the confirmation of another approved medical measurement procedure.
- Should the unit be used to control the liquid volume flow, the user has to analyze the risk of the application and, if necessary, take additional actions independent from the system to minimize this risk.

## 5.7 Electromagnetic Compatibility

Medical electronic devices are subjected to special precautions in terms of the electromagnetic compatibility (EMC). Because of this, the SonoTT™ Clamp-On Transducer and SonoTT™ Ultrasonic FlowComputer have to be installed and operated according to the following EMC information.

The SonoTT™ Ultrasonic FlowComputer and the SonoTT™ Clamp-On Transducer fulfill the requirements of the *IEC 60601-1-2 4<sup>th</sup> ed.*. Compliance with these requirements can only be guaranteed if the device is used according to the operating instructions.

This device fulfills the requirements of Class A according to CISPR 11 and should not be operated in home care, industrial areas, or environments with strong electromagnetic fields or disturbances.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the SonoTT™ Clamp-On Transducer or SonoTT™ Ultrasonic FlowComputer including the cables. Otherwise, degradation of the measurement performance and stability of this equipment could result.

If an electric cautery comes in contact with the device or the SonoTT™ Ultrasonic FlowComputer, it might cause a dysfunction which can lead to incorrect measurement values or other problems. This includes use of HF surgery at the patient or crossing of lines. Increased emission or decreased electromagnetic noise immunity can affect the measurements.

They can be caused by the use of SonoTT™ Clamp-On Transducers or any devices and/or accessories that are part of the measurement system but which are not mentioned in these operating instructions. Interference of ultrasonic flow measurements by electromagnetic fields could be possibly identified by compromised measured data, which are not related to the real flow. Please also note that the SonoTT™ Clamp-On Transducer is not galvanically isolated.

The SonoTT™ Clamp-On Transducer performs without producing either unintended or excessive ultrasound output nor unintended or excessive surface temperatures according to *IEC 60601-1-2 4<sup>th</sup> ed.* Electromagnetic fields and electrostatic discharge cannot lead to unintended changes of the electronics causing hazardous ultrasonic output. Nevertheless, it is important to ensure that the SonoTT™ Clamp-On Transducer and the SonoTT™ Ultrasonic FlowComputer are not placed near any disturbance source that is not compliant with the applicable standards, since this could influence and negatively affect the measurement.

## 5.8 Ultrasonic Safety Regarding Acoustic Output

When in use, the connected SonoTT™ Clamp-On Transducers produce very low-power ultrasonic signals at frequencies between 2.25 and 4 MHz. The emitted acoustic signals are declared in the technical specifications according to *IEC 61157 2<sup>nd</sup> ed.* and *IEC 60601-2-37 2<sup>nd</sup> ed.*

If used properly, the acoustic signals produced by the SonoTT™ Ultrasonic FlowComputer and SonoTT™ Clamp-On Transducer are, according to current knowledge, harmless for people and the environment.

## 5.9 Biocompatibility

The SonoTT™ Clamp-On Transducer is made from conventional materials which, with appropriate use, are harmless for the patient and the user according to current knowledge.

## 5.10 Maintenance and Service

The service for the SonoTT™ Clamp-On Transducer may be carried out by em-tec GmbH only. If these instructions are not followed, em-tec GmbH shall accept no liability for the component and the warranty will be void.

In order to maintain the accuracy of the SonoTT™ Clamp-On Transducer, it is recommended to regularly carry out re-calibrations. The time between calibrations should not exceed 24 months.

If you experience any trouble with the measurement despite following the operating instructions, or if your SonoTT™ Clamp-On Transducer is damaged in any way, please contact our service department.

Please make a note of the serial number of the SonoTT™ Clamp on Transducer and the evaluation device before you contact our staff.

If you need to return the SonoTT™ Ultrasonic FlowComputer or the SonoTT™ Clamp-On Transducer for servicing, please follow the steps below:

1. Contact our service department at:  
em-tec GmbH  
Service Department  
Am Graben 6-8  
86923 Finning  
Germany  
em-tec-service@psgdover.com
2. Our service department will send you a RMA form.
3. Fill out the form and include it in your shipment.

## 6 Installation, Use, and Disassembly of the SonoTT™ Clamp-On Transducer

### 6.1 General Safety Information

- If the SonoTT™ Clamp-On Transducer is not used as intended, the user may be exposed to risks that were not taken into account during its development.
- If the SonoTT™ Clamp-On Transducer is clamped onto a tubing circuit that is part of a life-supporting device used to maintain a patient's bodily functions, to support the lives of critically-ill patients, or during surgical interventions, these devices must be used and monitored by qualified clinical professionals.
- Decisions regarding therapeutic methods must not be solely based on the measured flow data.
- The values issued by the flow measurement system including the SonoTT™ Clamp-On Transducer have to always be verified with other methods before any decisions regarding treatment are made.
- For the installation, all ambient conditions must meet the prescribed specifications (see chapter 10 “Technical Specifications of the SonoTT™ Clamp-On Transducer”).

### 6.2 Connecting the SonoTT™ Clamp-On Transducer to the SonoTT™ Ultrasonic FlowComputer

To connect the SonoTT™ Clamp-On Transducer to the SonoTT™ Ultrasonic FlowComputer, carefully attach the 15-pin high density D-Sub socket to the screw connection at the back of the SonoTT™ Ultrasonic FlowComputer. Once this is done, the sensor can be attached to the tube.

- The connector of the SonoTT™ Clamp-On Transducer must be dry when plugged in.
- Check the plug and the pins for damage before connecting the sensor.
- Do not use the SonoTT™ Clamp-On Transducer if any part of it is damaged.
- Check and make sure that the connector used fits and is actually compatible with the measurement system.

## 6.3 Attaching the SonoTT™ Clamp-On Transducer to the Tube

EN

- Before you clamp on the sensor, make sure that the tube is not deformed, damaged, or dirty.
- Air in the tube can lead to errors in the measurement or interrupt the ultrasound coupling.
- Ensure that there is no debris or dirt in the flow channel. If necessary, clean the flow channel before attaching the SonoTT™ Clamp-On Transducer to the tube.
- While installing the SonoTT™ Clamp-On Transducer, be sure not to impair the function of the tubing system.
- Ensure that the tube size and material used corresponds with the size and material the SonoTT™ Clamp-On Transducer was calibrated and adjusted for.
- First, make sure that the arrow on the lid is aligned with the flow direction (positive flow display).

**Note:** If the medium within the tube moves in the opposite direction of the arrow on the lid, the flow value is displayed as negative value (e.g. -250 ml/min instead of 250 ml/min). The accuracy and the flow range stated in chapter 10 "Technical Specifications of the SonoTT™ Clamp-On Transducer" remains the same for both positive and negative values.

To clamp on the SonoTT™ Clamp-On Transducer, follow the steps illustrated below:


Figure 3: Clamping on the Sensor

1. Open the spring-loaded lock and swing back the lid.
2. Insert the tube into the channel.
3. Close the SonoTT™ Clamp-On Transducer lid and secure it with the help of the spring-loaded lock.

- Do not use excessive force when closing the sensor lid; while the tube should be fixed within the channel, it must not be crammed into it as this might affect the measurement values or damage the sensor.
- Ensure that the lid is fully closed as to avoid fluctuating coupling and flow values.
- To avoid flow turbulences and associated measurement inaccuracies, it is necessary for the tube to be straight in the area near the SonoTT™ Clamp-On Transducer.
- Ideally, the tube should be straight for a distance of at least 10 x the inner diameter (ID) of the tube on either side of the sensor.


Figure 4: Ideal Inlet Section

For more information regarding the setup of the SonoTT™ Clamp-On Transducer and what to keep in mind to maintain and optimize the accuracy of your flow measurement, refer to our TechNote "Optimizing Accuracy". The TechNote is available upon request, or for download on our website.

## 6.4 Use of the SonoTT™ Clamp-On Transducer

For more details on how to carry out the flow measurement, please refer to the operating instructions of the evaluation device; in this case, the SonoTT™ Ultrasonic FlowComputer.

EN

Make sure to follow the steps below and to heed the instructions included in them in order to ensure the highest possible accuracy of measurement values as well as the highest possible safety for the patient.

In addition, always ensure that

- the medium temperature is the temperature the sensor was calibrated and adjusted for.
- there is a zero flow adjustment performed before each measurement.
- the sensor is not permanently placed on the patient as this might lead to skin irritations.

### Before starting the flow measurement:

1. Make sure that only tubes specified for the SonoTT™ Clamp-On Transducer are used as this is the tube type the SonoTT™ Clamp-On Transducer is adjusted and calibrated for.
2. Before each measurement, the cables, connectors and pins must be checked for breaks or damages as this could result in wrong measurements. If there are any damages, do not use the SonoTT™ Clamp-On Transducer.
3. Ensure that the temperature of the medium is the one the SonoTT™ Clamp-On Transducer was adjusted and calibrated for.
4. Ensure that no air bubbles are in the tube during the zero flow adjustment. Allow sufficient time for the system to adapt to ambient conditions, then zero the flow. Only then is the system ready for measurement.
5. The coupling must be checked before each measurement.
6. **A zero flow adjustment must be carried out before each measurement** to avoid possible offsets from the measured values.

**Before the zero flow adjustment**, allow the system sufficient time to adapt to ambient conditions.

#### During the zero flow adjustment,

- The sensor must be clamped on the correct position .
- The tube must be filled with liquid, but the liquid must not move.
- There must be no air bubbles inside the tube.

### During the flow measurement:

- Ensure that the coupling is within the operating range of 60 % - 99%. If the coupling falls below 50 %, the SonoTT™ Ultrasonic FlowComputer issues a warning as the low coupling might impair the measurement.
- It is important that any unwanted interruption of the extracorporeal patient circulation due to damage or kinking of the tubing is avoided at all times.
- The stated measurement accuracy of the SonoTT™ Clamp-On Transducer is only valid if the specified tubes are used and if the temperature of the medium is the one the sensor was adjusted and calibrated for.
- Any changes in the medium itself or its temperature can cause errors or anomalies in the measurement.
- It must be ensured that the correct calibration parameters are used.

While it is possible to carry out measurements outside of the specified flow range, em-tec GmbH cannot guarantee that resulting flow values will be within the given accuracy range.

The parameters of the connected SonoTT™ Clamp-On Transducer and the selected calibration table are displayed in the **Settings-Menu** of the SonoTT™ Ultrasonic FlowComputer as shown in the following example:

Calibration Table	Transducer Serial No.	Tube Size	Tube material	Medium	Medium Temperature	Maximum Measured Value
0	SN 12345	3/8" x 3/32"	PVC	Blood (based on water)	37 °C (98.6 °F)	10000 ml/min

Please ensure that the tube size and material used corresponds to that shown on the display of the SonoTT™ Ultrasonic FlowComputer.

## 6.5 Disassembly and Storage of the SonoTT™ Clamp-On Transducer

In order to not impair the measurement, it is best to disassemble the SonoTT™ Clamp-On Transducer once the measurement process has been completed.

To disconnect the SonoTT™ Clamp-On Transducer, follow these steps:

1. Open the spring-loaded lock and swing back the lid.
2. Remove the tube from the channel.
3. Close the lid and carefully remove the 15-pin, high density D-Sub socket from the screw connection at the back of the SonoTT™ Ultrasonic FlowComputer.

Once you have disconnected the SonoTT™ Clamp-On Transducer, check for any residue of the medium or any other kind of dirt or grime and clean the sensor if necessary. Store the sensor according to storage specifications. Make sure to store it in a place where it cannot fall off or be dropped and ensure that the plug is dry.

## 6.6 Cleaning and Disinfection of the SonoTT™ Clamp-On Transducer

Make sure no detergents, disinfectants or other liquids can leak into the plug of the SonoTT™ Clamp-On Transducer.

For the cleaning of the device, we recommend the use of common household detergents. In order to avoid stains and dirt from drying out, grime residue should be removed after each application.

To clean the SonoTT™ Clamp-On Transducer, use a lint-free cloth. Be sure to remove any dirt in the clamping area. Additionally, please follow the legal regulations valid in your country and the hygiene regulations for your doctor's office or hospital.

- The SonoTT™ Clamp-On Transducer may not be submerged into cleaning or disinfecting solution.
- The SonoTT™ Clamp-On Transducer is not suited for cleaning processes using machines. Sterilization processes, especially steam sterilization or autoclaving, may not be used.
- The SonoTT™ Ultrasonic FlowComputer and the SonoTT™ Clamp-On Transducer are used in non-sterile areas in operating rooms. Due to this, cleaning and disinfection after every use as well as before the first use is essential.

The following surface disinfectant is recommended:

Name	Manufacturer	Contact
Bacillo® AF	Hartmann	www.hartmann.de

## 6.7 Changing the SonoTT™ Clamp-On Transducer

If the SonoTT™ Clamp-On Transducer was exchanged before switching on the SonoTT™ Ultrasonic FlowComputer, then the new calibration parameters will be shown in the Settings Menu of the SonoTT™ Ultrasonic FlowComputer after the next start-up sequence.

If the SonoTT™ Clamp-On Transducer is exchanged during the operation of the SonoTT™ Ultrasonic FlowComputer, the alarm is activated with an error message. Once the new SonoTT™ Clamp-On Transducer is connected, the alarm is reset and the calibration data can be seen in the Settings Menu.

- Make sure that the SonoTT™ Clamp-On Transducer is attached at a different place on the tube each time to avoid deformation.
- Ensure that the calibration table used is the one needed for your respective application.

## 7 Adjustment and Calibration of the SonoTT™ Clamp-On Transducer

As the function of the SonoTT™ Clamp-On Transducer is based on an acoustic measurement principle, parameters such as the tubing size and material, or the medium as well as its temperature can impact the ultrasonic signal and thus the transit time difference between them. To take these factors into account and to ensure the highest possible accuracy regarding the measured values, each sensor is individually adjusted and calibrated prior to shipment.

The relevant parameters for the adjustment and calibration are:

- tubing material
- medium
- medium temperature
- flow range

### 7.1 Adjustment

An act of modifying something, in this case to fit a certain standard, the sensor adjustment represents an actual intervention in the measurement system. To ensure that the parameters stated above are factored into the measurement and consequently into the flow values displayed by the SonoTT™ Ultrasonic FlowComputer, parameters such as the amplitude of the ultrasonic signal or the length of the ultrasonic burst are changed to fit the tubing and the medium used during the application. .

As a result, the adjustment is an action that permanently alters the device as the setting of the burst signal (scan window and amplitude) and the sequential balance line are changed.

Since these parameters have a deep impact on the sensor's performance, and since adjusting the sensor represents an intervention of the measurement system, adjustments may only be carried out by em-tec GmbH.


Figure 5: Ultrasonic Signal

### 7.2 Calibration

The term 'calibration', on the other hand, refers to the act of checking the accuracy of a measurement device by comparing its measurement values with those of a device with a known accuracy. This comparison is subject to a specified measurement setup and specified measurement conditions. During the calibration process, neither the devices involved nor any of their parameters may be changed, manipulated, or altered in any way.

Ideally, the SonoTT™ Clamp-On Transducers should be re-calibrated every 24 months. This can be carried out by the customer on-site or by em-tec GmbH.

For more information regarding this process, please contact em-tec GmbH:

## 8 Troubleshooting

If there are any issues with the SonoTT™ Clamp-On Transducer, try the following suggestions. If the problem cannot be solved this way, please contact your local distributor or em-tec GmbH directly.

Problem	Possible Cause	Action
Sensor lid does not close.	<ul style="list-style-type: none"> <li>• Sensor might be damaged.</li> <li>• Tube might be too large.</li> </ul>	<ol style="list-style-type: none"> <li>1. Check if the sensor is damaged.</li> <li>2. Check if the tube size is the one indicated for your sensor.</li> <li>3. If necessary, exchange the tube for the right one.</li> </ol> <p>If the problem persists, return the sensor for servicing*.</p>
Sensor does not stay in place.	<ul style="list-style-type: none"> <li>• Sensor might be damaged.</li> <li>• Tube might be too small.</li> </ul>	<ol style="list-style-type: none"> <li>1. Check if the sensor is damaged.</li> <li>2. Check if the tube size is the one indicated for your sensor.</li> <li>3. If necessary, exchange the tube for the right one.</li> </ol> <p>If the problem persists, return the sensor for servicing*.</p>
Sensor emits noise.	There might be electromagnetic fields influencing the sensor.	<ol style="list-style-type: none"> <li>1. Check if there are any electromagnetic fields or disturbance sources that are not compliant with the relevant standards and that could influence the sensor. Also check if the sensor or its cable came into contact with an electric cautery.</li> <li>2. Remove the source of the electromagnetic field or carry the measurement out somewhere else.</li> </ol> <p>If the problem persists, return the sensor for servicing*.</p>
Sensor is not recognized.	The cable is damaged.	<ol style="list-style-type: none"> <li>1. Check if the sensor is connected properly to the device and if the sensor is damaged.</li> <li>2. Check if the cable is damaged.</li> <li>3. If the cable is damaged, return the sensor and cable for servicing.</li> </ol> <p>If the problem persists, return the sensor for servicing*.</p>
There is no measurement possible.	The plug is damaged (e.g. broken PIN).	<ol style="list-style-type: none"> <li>1. Check if the sensor is connected properly to the device and if the sensor is damaged.</li> <li>2. Check if the plug is damaged and if there are any broken PINs</li> <li>3. If the plug is damaged, return the device for servicing.</li> </ol> <p>If the problem persists, return the sensor for servicing*.</p>

Problem	Possible Cause	Action
<p>Sensor is not recognized.</p> <p>There is no measurement possible.</p>	<p>Storage might have been damaged or deleted by electromagnetic fields.</p>	<ol style="list-style-type: none"> <li>1. Check if the sensor is connected properly to the device and if the sensor is damaged.</li> <li>2. Check if there are any electromagnetic fields or disturbance sources that are not compliant with the relevant standards and that could influence the sensor. Also check if the sensor or its cable came into contact with an electric cautery.</li> <li>3. Remove the source of the electromagnetic field or carry the measurement out somewhere else.</li> </ol> <p>If the problem persists, return the sensor to em-tec GmbH for recalibration*.</p>
	<p>The connected sensor is the wrong sensor type.</p>	<ol style="list-style-type: none"> <li>1. Check if the correct sensor is connected and if the connected sensor is not damaged.</li> <li>2. If the sensor is correct, disconnect it.</li> <li>3. Reconnect the sensor again.</li> </ol> <p>If the problem persists, return the device and sensor for servicing,</p>
<p>Acoustic coupling is too low.</p>	<p>Sensor lid might not be closed properly.</p>	<ol style="list-style-type: none"> <li>1. Check if the sensor lid is closed properly.</li> <li>2. Close the lid.</li> </ol> <p>If the problem persists, return the sensor for servicing*.</p>
	<p>The tube size and/or material is not the one indicated for the sensor.</p>	<ol style="list-style-type: none"> <li>1. Check if the tube size and material are the ones the sensor was adjusted and calibrated for.</li> <li>2. If necessary, exchange the tube for the right one.</li> </ol> <p>If the problem persists, return the sensor for servicing*.</p>
	<p>There is dirt in the flow channel.</p>	<ol style="list-style-type: none"> <li>1. Check if there is dirt in the flow channel.</li> <li>2. If necessary, clean the flow channel.</li> </ol> <p>If the problem persists, return the sensor for servicing*.</p>
	<p>Sensor might be damaged.</p>	<ol style="list-style-type: none"> <li>1. Check if the sensor or the sensor lid is damaged (e.g. broken hinge).</li> <li>2. If the sensor is damaged, stop the measurement and return the sensor for servicing.</li> </ol> <p>If the problem persists, return the sensor for servicing*.</p>
	<p>There are air bubbles in the tube.</p>	<ol style="list-style-type: none"> <li>1. Check if there are air bubbles in the tube.</li> <li>2. Remove any air within the tube.</li> <li>3. Repeat the flow measurement following the steps described in chapter 6.4.</li> </ol> <p>If the problem persists, return the sensor for servicing*.</p>

Problem	Possible Cause	Action
Measurement appears to be false.	Sensor might be damaged.	<ol style="list-style-type: none"> <li>1. Check if the sensor or the sensor lid (e.g. broken hinge) is damaged.</li> <li>2. If the sensor is damaged, stop the measurement and return the sensor for servicing.</li> </ol> If the problem persists, return the sensor for servicing*.
	Sensor lid might not be closed properly.	<ol style="list-style-type: none"> <li>1. Check if the sensor lid is closed properly.</li> <li>2. Close the lid.</li> </ol> If the problem persists, return the sensor for servicing*.
	The application parameters differ from the ones the sensor was adjusted and calibrated for (i.e. tube size or material, medium).	<ol style="list-style-type: none"> <li>1. Check if the application parameters are the ones the sensor was adjusted and calibrated for.</li> <li>2. If necessary, exchange the tube for the right one.</li> </ol> If the problem persists, return the sensor for servicing*.
	There is debris or dirt in the flow channel.	<ol style="list-style-type: none"> <li>1. Check if there is debris or dirt in the flow channel.</li> <li>2. If necessary, clean the flow channel.</li> </ol> If the problem persists, return the sensor for servicing*.
	There might be electromagnetic fields influencing the sensor.	<ol style="list-style-type: none"> <li>1. Check if the sensor is connected properly to the device and if the sensor is damaged.</li> <li>2. Check if there are any electromagnetic fields or disturbance sources that are not compliant with standards that could influence the sensor or if the sensor or its cable came into contact with an electric cautery.</li> <li>3. Remove the source of the electromagnetic field or carry the measurement out somewhere else.</li> </ol> If the problem persists, return the sensor to em-tec GmbH for recalibration*.
	Wrong calibration table selected.	<ol style="list-style-type: none"> <li>1. Check if the selected calibration table is the one appropriate for your application.</li> <li>2. If necessary, select the right calibration table.</li> </ol> If the problem persists, return the sensor for servicing*.

Problem	Possible Cause	Action
Measurement appears to be false.	The plug is damaged (e.g. broken pin).	<ol style="list-style-type: none"> <li>1. Check if the sensor is connected properly to the device and if the sensor is damaged.</li> <li>2. Check if the plug is damaged and if there are any broken PINs.</li> <li>3. If the plug is damaged, return the device for servicing.</li> </ol> <p>If the problem persists, return the sensor for servicing*.</p>

For more troubleshooting, please refer to the operating instructions of the SonoTT™ Ultrasonic FlowComputer.

\*For more information concerning the return and servicing of our products, please refer to chapter 5.9 "Maintenance and Service".

## 9 Environmental Protection and Disposal

### Packaging

Packaging materials are made from environmentally friendly materials. The packaging materials will be disposed of by em-tec GmbH upon request.

### Disposal

The SonoTT™ Clamp-On Transducer and its accessories must be disposed of in accordance with the applicable national provisions for electronic components. In accordance with the requirements of EU Directive **2012/65/EC** Waste Electrical and Electronic Equipment (WEEE), our customers in the EU are entitled to return all waste deriving from the products to us – in clean and disinfected condition. The em-tec GmbH WEEE registration number is: **DE 99135207**.

Upon receipt, we repair or dispose of these components properly. For our address please see the very beginning of these operating instructions. For the best utilization of raw materials, the product and its components and accessories should not be disposed of together with household waste. All parts must be collected separately from household waste and disposed of in an environmentally responsible way in accordance to local regulations.

- Before disposal, decontaminate all parts according to the procedure applicable in the clinic.
- If you have questions about disposal, please contact em-tec GmbH's service department.
- Waste may only be brought to the appropriate recycling facility if there is no risk of potential infection from electrical and electronic waste.

## 10 Technical Specifications of the SonoTT™ Clamp-On Transducer

Size (HxWxD), Weight (incl. cable & plug)	25 x 33 x 45 mm, 137 g (IDs 13479; 13481; 13475; 13478; 13480; 13476) 27 x 38 x 51 mm, 153 g (ID 13477)
Housing and Lid Material	Epoxy Resin, Aluminum, Brass
Cable length	2.9 m ±5 cm
Connector plug	15-pin high density D-Sub
Safety Class	Type BF applied part (in combination with the SonoTT™ Ultrasonic FlowComputer)
IP-Code	without connector: IP67 (protected against splashing water)
Recommended tube type	flexible, non-reinforced tubing; e.g. silicone, PVC
Acoustic output signal	Frequency 2,25 MHz, $p_r < 1 \text{ MPa}$ , $I_{\text{spta}} < 20 \text{ mW/cm}^2$ , $I_{\text{ob}} < 100 \text{ mW/cm}^2$ , $\Delta T < 0.3 \text{ }^\circ\text{C}$ , $TI < 1.0$ , $MI < 1.0$
Possible adjustment temperature range	4 °C to 41 °C
Expected product life:	10 years

### Accuracy<sup>1</sup> in Combination with the SonoTT™ Ultrasonic FlowComputer:

Zero flow deviation: ±15 ml/min

### Flow Accuracy<sup>2</sup>

0 < flow < Qmin	not defined
Qmin to 1000 ml/min	±100 ml/min
> 1000 ml/min	±7 % of the value ±30 ml/min

<sup>1</sup>The accuracy stated above can only be guaranteed if the parameters throughout the measurement are the same as the ones the sensor was adjusted and calibrated for (e.g. medium type, medium temperature, tube size, tube material).

### <sup>2</sup> Please note:

- Qmin refers to the minimum flow value for which the accuracy is specified.
- Qmax refers to the maximum flow value for which the accuracy is specified.
- The given values describe the positive flow range, but are the same for negative flow values.

### Transport and Storage Conditions:

Atmospheric Pressure	70 kPa to 106 kPa
Temperature Range	-20 °C to 55 °C (-4 °F to 131 °F)
Relative Humidity	10 % to 96 % (non-condensing)

### Ambient Conditions during Operation:

Atmospheric Pressure	70 kPa to 106 kPa
Temperature Range	10 °C to 40 °C (50 °F to 104 °F)
Relative Humidity	10 % to 96 % (non-condensing)

## Range of SonoTT™ Clamp-On Transducers and Their Flow Measurement Range

ID	Type	Qmin <sup>3</sup>	Qmax <sup>4</sup>	Tube Size (Inner diameter (ID) x Wall Thickness (WT))
13482	CT 5,7 mm	100 ml/min	±1000 ml/min	5.7 mm outer diameter
13479	CT 6,8 mm	100 ml/min	±6000 ml/min	6,8 mm outer diameter
13481	CT 3/16" x 1/16"	100 ml/min	±6000 ml/min	3/16" x 1/16"
13478	CT 1/4" x 3/32"	120 ml/min	±8 000 ml/min	1/4" x 3/32"
13475	CT 1/4" x 1/16"	120 ml/min	±8 000 ml/min	1/4" x 1/16"
13476	CT 3/8" x 3/32"	150 ml/min	±10 000 ml/min	3/8" x 3/32"
13480	CT 3/8" x 1/16"	150 ml/min	±10 000 ml/min	3/8" x 1/16"
13477	CT 1/2" x 3/32"	300 ml/min	±20 000 ml/min	1/2" x 3/32"

<sup>3,4</sup> **Please note:** While a measurement is possible for the flow ranges stated above, the accuracy is only defined for the flow range between the Qmin and the Qmax as stated under "Accuracy in Combination with the SonoTT™ Ultrasonic FlowComputer".

**Please note:** If the medium within the tube moves in the opposite direction of the arrow on the lid, the flow value is displayed as negative value (e.g. -250 ml/min instead of 250 ml/min). The accuracy and the flow range stated above is the same for positive as well as negative flow values.

Please specify calibration parameters at time of ordering, i.e. media, temperature and tubing type.

## About em-tec GmbH

em-tec has been a specialist for flow measurement systems in the medical and bioprocessing technology sector for over 30 years. The company's core competence is the non-invasive flow measurement using the ultrasonic transit-time method, that is used for applications in extracorporeal circulation systems of life-sustaining systems as well as in biopharma applications that use flexible tubes. Headquartered in Finning, Germany, em-tec is part of PSG®, a Dover company.

For more information about em-tec, please visit [em-tec.de](http://em-tec.de). For more information about PSG®, please visit [psgdover.com](http://psgdover.com).

## About Dover

Dover is a diversified global manufacturer and solutions provider with annual revenue of approximately \$7 billion. We deliver innovative equipment and components, consumable supplies, aftermarket parts, software and digital solutions, and support services through five operating segments: Engineered Products, Fueling Solutions, Imaging & Identification, Pumps & Process Solutions and Refrigeration & Food Equipment. Dover combines global scale with operational agility to lead the markets we serve. Recognized for our entrepreneurial approach for over 60 years, our team of over 23,000 employees takes an ownership mindset, collaborating with customers to redefine what's possible. Headquartered in Downers Grove, Illinois, Dover trades on the New York Stock Exchange under "DOV." Additional information is available at [dovercorporation.com](http://dovercorporation.com).

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Für weitere Informationen zu diesem Verfahren wenden Sie sich bitte an em-tec GmbH.
- Wenn Sie Fragen zur Entsorgung haben, wenden Sie sich bitte an die Service-Abteilung der em-tec GmbH.
## Über em-tec GmbH